


Advancing the Clinical Detection of Hemostatic Abnormalities

CHRONO-LOG Model 700

Whole Blood/Optical Lumi Aggregometer

- **Easy to Use with Rapid Results**
- **Whole Blood and Optical Aggregation Modes**
- **Disposable and Reusable Impedance Electrodes**
- **Simultaneous ATP Secretion Studies**
- **Dual and Four Channel Systems**



A Proven Track Record of Accurate and Reliable Results

The Answer to Your Platelet Function Testing Needs

Rapid and Simple Screening

Chrono-log Whole Blood/Optical Lumi-Aggregation Systems improve detection of abnormal platelet function by simultaneously measuring platelet aggregation and dense granule secretion in a physiologic whole blood environment or traditional plasma sample. In addition to optical aggregation, the turbidometric feature permits platelet function testing for the Ristocetin Cofactor Assay and Sticky Platelet Syndrome tests.

Efficient and economic, the simplicity and heightened sensitivity of whole blood testing in the Clinical Laboratory permits the establishment of a therapeutic monitoring program:

- Monitor aspirin therapy to determine response or non-response with Collagen-induced aggregation in whole blood.³⁸⁻⁴¹
- Monitor Plavix therapy with ADP-induced aggregation in whole blood.^{42,43}
- Monitor administration of DDAVP by Ristocetin-induced aggregation in whole blood.^{2,4,8}
- Measure effect of GPIIb/IIIa receptor blockage and other anti-platelet drugs.³¹

Identify Patients at Risk of Thrombosis or Bleeding... Select Patients for Anti-Platelet Therapy... Monitor Dosage and Response

Patients with a previous history of thrombosis show a high prevalence of platelet hyperactivity. Where the goal of therapy is inhibiting platelet aggregation, Chrono-log Whole Blood Aggregometers can routinely identify patients at risk for thrombosis or bleeding and can verify effective treatment.³⁷ For example:

Aspirin (ASA) is used therapeutically to inhibit thromboxane production of circulating platelets. Measuring platelet aggregation in whole blood with low and high Collagen concentrations (1 and 5 ug/mL) is a cost-effective, quick screen for monitoring effective dosage levels of ASA. While low Collagen concentration at 1ug/mL can best detect the aggregation defect caused by ASA,³⁶ the high Collagen concentration of 5ug/mL bypasses the ASA effect, serving as an auto-control and allowing the level of response to be determined. Poor aspirin response is defined as $\geq 50\%$ of platelets aggregating when the 1ug/mL is compared to the 5ug/mL test.³⁸⁻⁴¹ Simultaneous measurement of thromboxane-dependent ATP release with Luminescence shows reduced ATP release with both Collagen concentrations.¹⁸

Plavix (Clopidogrel) inhibits the binding of ADP to its receptor and the subsequent ADP-mediated activation of the glycoprotein GPIIb/IIIa complex, inhibiting platelet aggregation. When testing in whole blood with 20 uM ADP, a 57-60% mean dose dependent inhibition of platelet aggregation was seen one day after a loading dose of 600 mg Clopidogrel was administered.⁴³ Following ten days of treatment with 75 mg Clopidogrel, whole blood aggregation with 2.5 and 5 uM ADP was absent while tests with ADP at 10 and 20 uM showed "markedly impaired aggregation." Optical aggregometry with PRP and ATP secretion remained normal; therefore, Whole Blood Aggregometry with ADP is an effective screen to determine a patient's response or non-response to Plavix.⁴²

The GPIIb/IIIa receptor antagonist abciximab (ReoPro™) inhibits platelet function by preventing the binding of fibrinogen and von Willebrand factor to activated platelets. Whole Blood Impedance aggregation with 5ug/mL Collagen more closely parallels receptor blockage than turbidometric aggregation with 5 or 20 uM ADP and is a rapid, simple and accurate method for measuring the effects of abciximab.³²

Dipyridamole (Persantin), a coronary vasodilator, anti-thrombotic drug, significantly inhibits aggregation in whole blood when testing with 3.5 uM ADP or 5 ug/mL Collagen; this inhibition is not seen in platelet-rich plasma. Interaction between platelets and red blood cells plays an important role in the anti-platelet activity of Dipyridamole.

In a study noted above, it was concluded that: "The WB impedance assay should be widely used in the clinical practice to identify subjects who are resistant to aspirin and clopidogrel."^{37, 42}

CHRONO-LOG Model 700

Whole Blood/Optical Lumi Aggregometer

Measures platelet function on patient samples using electrical impedance in whole blood or optical density in plasma while simultaneously measuring ATP release by the luminescence method.

Comprehensive Diagnostic Capability

Platelet disorders shown below can be detected using Whole Blood/Optical Aggregometry and/or Luminescence:

von Willebrand Disease	Extracorporeal circulation	Deficiency of the enzyme cyclo-oxygenase
Glanzmann's Thrombasthenia	Heparin-Induced Thrombocytopenia	Deficiency of the enzyme thromboxane synthetase
Storage Pool Deficiency (SPD)	Sticky Platelet Syndrome	Thrombocytopenia with absent radii (TAR syndrome)
Thrombopathia or Thrombocytopenia	Bernard-Soulier Syndrome (BSS)	Risk of thrombosis
May-Heggelin Anomaly	Gray Platelet Syndrome	Membrane receptor site defects
Defective early responses (primary defects)	Secretion defects	
Non-steroidal, Anti-inflammatory drugs	Hyperaggregability	

Patients with histories suggesting hemorrhagic disorders are often investigated for coagulation defects rather than platelet dysfunction,⁷ because optical aggregation studies are time-consuming and ambiguous. NOW, with less than 5mL of blood and in less time than it takes to prepare plasma for tests of coagulation,^{5,6} you can detect a wide variety of platelet disorders.

For detection of vW disease, Ristocetin-induced platelet aggregation (R.I.P.A.) by the impedance method in whole blood, is a highly sensitive and time efficient method of screening persons at risk. While patients with von Willebrand disease exhibit the same abnormalities in whole blood as in platelet rich plasma, an extended lag phase of >70 seconds in whole blood allows for a clearer separation of vW patients from normals.^{7,8} In addition,

screening for Type 2B and Platelet-Type von Willebrand is quickly performed with low concentration Ristocetin.^{4,4}

As impedance aggregometry in whole blood is a more adequate tool for the detection of platelet hyperaggregability than the optical method in PRP,¹² more patients can be identified as at risk for thromboembolic complications. This greater sensitivity results from performing tests in the presence of red cells and leukocytes.¹³⁻¹⁵

The electronic impedance technique can measure aggregation when optical methods cannot — in hemolysed, icteric or lipemic samples where the sample turbidity interferes with measurement, or in "giant" platelet syndromes where centrifugation results in the loss of the platelets to be studied.^{28,29} Due to the small blood draw required for Whole Blood testing, [5 mL vs. 20 mL],

pediatric patients can be tested for platelet disorders. Specific platelet defects have been diagnosed with platelet counts as low as 70,000/mm³ with pediatric patients.³

Direct measurement of ATP secretion via Luminescence in whole blood or PRP provides unequivocal evidence of normal dense granule release and is a quantitative measure of platelet activation.⁹ In contrast, repetitive dose response optical aggregation tests in PRP failed to diagnose 67% of patients with Storage Pool Deficiency and a prolonged bleeding time.¹⁰ Simultaneous measurements of the release reaction provides further insight into the mechanisms of platelet response.¹¹

The CHRONO-LOG Model 700

Easy to Use... Rapid Results... Fast Turnaround Time

- Test Four Samples at One Time with Two Dual Channel Modules
- Two aggregation modes—Whole Blood (Impedance) and Optical (Turbidometric)
- Disposable and Reusable Impedance Electrodes
- ATP Secretion Studies (Luminescence)
- Push Button Controls
- LCD Display with System Error Monitoring for each Channel
- Built in Computer Interface
- AGGRO/Link8 and vW Cofactor Software Packages Included



CHRONO-LOG Model 700

Whole Blood/Optical Lumi Aggregometer

Efficient and Economic

Chrono-log Whole Blood/Optical Lumi-Aggregation Systems allow the study of platelet function in whole blood in the presence of other cells, before the decay of labile modulators. A micro-sample of whole blood is all that is required for measurement of aggregation and secretion, making the Chrono-log System particularly suitable for research with smaller laboratory animals.²⁰⁻²¹

The simultaneous measurement of aggregation with a secretion study simplifies, speeds and reduces the cost of your platelet function testing. In addition to saving time, less blood is needed, precious agonists are conserved and consumable costs reduced.

A Powerful and Flexible Research Tool (Three Systems in One)

• Optical Aggregation

Optical (turbidometric technique)

The Chrono-log Model 700 Aggregation Systems include automatic optical aggregation channels for measuring platelet or leukocyte aggregation in PRP or isolated samples, for measuring Ristocetin CoFactor Activity, and for measuring agglutination of latex particles or investigating shape change of suspended cells.


Dual-beam infrared light sources and photodiode detectors are used to determine the difference in light transmission between PRP and PPP, with the zero and full-scale baselines automatically set by a single pushbutton. A platelet count difference of only $50 \times 10^9/L$ is needed for a full-scale deflection. Laboratory personnel can check the calibration of the optical

circuitry and perform an auto-calibration, greatly reducing down time. This feature provides more reliable and accurate calculation, which is especially important in the calculation of Slope in the performance of the Ristocetin CoFactor Assay.

• Impedance Aggregation



Disposable Probe



Reusable Probe



• Luminescence

Sensitive Luminescence Technique

A sensitive photomultiplier tube in the Chrono-log Lumi-Aggregation Systems provides a voltage output proportional to the luminescence intensity. A wide variety of luminescent reactions (ATP release, total adenine nucleotides, superoxide generation and others) can be measured by a Chrono-log Lumi-Aggregation System.

Bio-Luminescence and Chemi-Luminescence Assays

- Adding CHRONO-LUME[®] Reagent (buffered firefly luciferin-luciferase) to the sample provides the ability to measure platelet ATP Release. In addition to detecting Storage Pool and Secretion defects in the Clinical Laboratory, Luminescence is also a rapid and sensitive method for detecting Heparin-Induced Thrombocytopenia (HIT).³³

- Addition of Luminol to a sample of Leukocytes provides the ability to measure Superoxides. Chemi-luminescence has been proven to be of value in the identification of neutrophil defects such as Chronic Granulomatous Disease and myeloperoxidase deficiency.¹⁷

White Blood Cells

[In the presence of leukocytes, platelet aggregation is amplified by the physiologic stimulus, platelet activating factor (PAF).]

- Chrono-log Whole Blood/Optical Lumi-Aggregometers can also be used for leukocyte function studies.^{16,34,35} Leukocyte aggregation and superoxide generation can be measured simultaneously in a single sample of diluted blood or prepared leukocytes.

- Aggregation has also been demonstrated when leukocytes were activated by N-formyl-methionyl-leucylphenylalanine (fMLP) and leukotriene B₄ (LTB₄).²³ These data add to the possibility that the onset and propagation of inflammation, hemostasis, thrombosis and atherosclerosis might be influenced by platelet-leukocyte interactions.²⁴⁻²⁶

CHRONO-LOG Model 700

Whole Blood/Optical Lumi Aggregometer

Specifications

Chrono-Log Patented Electronic Impedance Technique

Chrono-log Whole Blood/Optical Lumi-Aggregation Systems detect platelet aggregation by passing a minute amount of electrical current between two electrodes immersed in whole blood, diluted blood, platelet rich plasma or washed platelets and by measuring the impedance between the electrodes. During initial contact with the sample, a monolayer of platelets forms on the electrodes. When an agonist is added

there is gradual platelet aggregation at the electrodes, increasing the impedance.²⁷ The increase in impedance upon stimulation is directly proportional to the platelets aggregating onto the Electrode Probe Assembly.²²

The Electrode Probe Assemblies have two carefully spaced rigid precious metal wires on which aggregation occurs. The

Reusable probe (P/N 369) has a plug-in cable long enough to permit easy cleaning between samples by rinsing in water and saline. The Disposable probe is supplied in its own test cuvette and is disposed along with the sample, test cuvette and stir bar following completion of the test. P/N 315 Disposable Probes are for single use only and do not require cleaning.

The 700 Models include both the Reusable Probes (P/N 369) and Disposable Probes (P/N 315)

Test Channels:

(2) or (4) Channels (Two, 2-Channel Modules), with:

Impedance Aggregation - aggregation in a 1mL sample. Automatic baseline set at 0% and 20Ω gain fixed at 50%.

Optical Aggregation - aggregation in PRP, gel-filtered or washed platelet samples.

Luminescence - Photo-multiplier tube detects ATP Release. Nine gain settings - X 0.005 to X 2.

Front Panel Display and Controls:

LCD Display - 24-characters x 2-line Liquid Crystal Display, one per channel, displays:

- Heater Block - temperature in °C
- Luminescence gain
- Stirring speed in RPM
- Operating mode (Impedance or Optical)
- Warning messages

Power ON/OFF Switch

Set Baseline Pushbutton(s) - Sets aggregation baselines to 0%. Adds 20 ohms (± 0.2 ohms) of resistance to set the impedance gain fixed at 50% of scale. Luminescence baseline set at 0%.

Select Switch - select Luminescence Gain, Temperature or Stirring Speed

Set Switch - Set Luminescence Gain, Temperature or Stirring Speed.

Mode Switch - Set Impedance or Optical mode

PPP Selector Switch - Set to 1, tests referenced to Channel 1 PPP. Set to own channel (2,3,4), test referenced to PPP for that channel.

Calibration Switch - Key-activated calibration of optical circuits.

Heater Block - set between 35.0°C and 39.0°C in 0.1°C steps. Error detection prevents operation when outside $\pm 0.2^\circ\text{C}$.

Stirrer - 400 to 1200 RPM in 100-RPM steps with "Stirrer Stopped" position. Error detection prevents operation if not within ± 10 RPM.

Sample Volumes

Whole Blood Lumi-Aggregation - typically 450μL whole blood plus 450μL of irrigation saline and 100μL CHRONO-LUME® Reagent.

Whole Blood Aggregation - typically 500μL whole blood plus 500μL irrigation saline

Optical Lumi-Aggregation - typically 450μL PRP plus 50μL CHRONO-LUME® Reagent; 225μL PRP with spacers

Optical Aggregation - typically 500μL PRP; 250μL with spacers

General Specifications (each module)

Power requirements - Switch selectable 115 or 230VAC ($\pm 10\%$), 50/60 Hz, 150 watts max.

Dimensions - 14" (36cm) wide, 8.5" (22 cm) high, 18" (46cm) deep

Weight - 40 lbs. (18kg)

Incubation Wells - Six (6) wells each channel @ $36.5^\circ \pm 1.0^\circ\text{C}$ when temperature set at 37°C .

Output Options:

Computer Interface - Digital Outputs - RS-232 and USB with AGGRO/LINK8 software.

Data Reduction System - (Included with Models 700-2DR and 700-4DR) - State of the Art Computer, Color Monitor and Color Printer. Software packages, Installed:

WINDOWS XP [subject to change]

AGGRO/LINK8 - real time color display of 4 channels of aggregation and ATP release... **total of eight traces**. Computes amplitude, slope, lag time and area under the curve. Reagent data, test values, and demographics stored in computer for later recall. [Requires *Windows 98 or better operating system*]

vW Cofactor Software - on-screen instructions, best-fit standard curve and CD calculated for 2 to 6 points. Allows rerun or deletion of serial point. Stores curves with lot numbers. Four samples can be run simultaneously. Percent of vW activity calculated and reported.

Chart Recorder - Analog output. Chrono-log Model 706 Single or Model 707 Dual Pen recommended. Other recorders must have (1) megohm min. input impedance, 100 mV range and 1 cm/minute chart speed.

Accessories and Supplies - Impedance

Electrode Probe Assembly [reusable] - P/N 369W (one for each channel)

Electrode Probe Assembly [disposable] - P/N 315-25, P/N 315-50, P/N 315-100 (includes test cuvette)

Cuvettes - 1mL, P/N 367

Stir Bars - Disposable: siliconized, P/N 370; Reusable: Magnetic, Teflon-Coated, P/N 368

Accessories and Supplies - Optical

Insert Assembly - P/N 366 (one for each channel)

Cuvettes - 500μL, P/N 312

Stir Bars - Disposable: siliconized, P/N 311; Reusable: Magnetic, Teflon-Coated, P/N 313

Spacers, P/N 365, for testing microvolume samples as low as 250μL

U.S. Patent No 6,004,818; European Patent No.

1004020 protects these instruments; other Patents Pending

Windows is a trademark of Microsoft Corporation

References

1. Riess H, Braun G, Brehm G and Hilkr E.: Critical evaluation of platelet aggregation in whole human blood. *Am J Clin Pathol*, Vol. 85, No. 1, (50-56) 1986.
2. Ibid
3. Bandi E, Thaiss H, Choi S, Kunzer W and Sutor AH: Platelet function test in citrated whole blood. Application and interpretation in pediatric patients. Auszug aus: Abstracts-Alexander-Schmidt-Gedachnis-Vorlesung 31. Jahrestag der Gesellschaft für Thrombose- und Hamostaseforschung (GTH) Freiburg i. Br. 18.-21. February 1987.
4. Sweeney JD, Hoening LA, Fitzpatrick JE: Whole blood aggregation in von Willebrand disease. *Am J Hematol*, 32, (190-193) 1989.
5. Ingerman-Wojenski CM and Silver MJ: A quick method for screening platelet dysfunctions using the whole blood lumi-aggregometer. *Thromb Haemost*. 51-2 (154-156) 1984.
6. Sweeney JD, Hoernig LA, Michnick A and Fitzpatrick JE: Whole Blood Aggregometry. Influence of sample collection and delay in study performance on test results. *Am J Clin Pathol*, Vol 92, No 5 (676-679) 1989.
7. Sweeney JD, Labuzetta JW, Bernstein ZR, Bielak KL and Fitzpatrick JE: Ristocetin-induced platelet aggregate formation and adherence to the probe of an impedance aggregometer. *Amer. Jnl. of Clin. Path.*, Vol. 93, No. 4, pp. 548-551, April 1990.
8. Sweeney et al, 1989.
9. Feinman RD, Derwiler TC and Ingerman-Wojenski C: The lumi-aggregometer as a research and clinical tool. *The Pits.: Physiology & Pharmacology*, (429-440) Copyright 1985 by Academic Press, Inc.
10. Nieuwenhuis HK, Akkerman J-WN and Sixma JJ: Patients with a prolonged bleeding time and normal aggregation tests may have Storage Pool Deficiency: studies on one hundred six patients. *BLOOD*. Vol. 70, No. 3, (620-623) September 1987.
11. Feinman et al
12. Abbate R, Boddi M, Frisco D and Gensini GF: Ability of whole blood aggregometer to detect platelet hyperaggregability. *Am J Clin Pathol*, Vol. 91, No. 2 (159-164) February 1989.
13. Ibid
14. Lacoste L, Latour JG, Theroux P and Waters DD: Hyperagregabilité plaquettaire et leucocytaire dans l'angine instable. *Med Sci (Suppl 1)*: 39A 1988.
15. Latour JG, Lacoste L, Theroux P, Leger-Gauthier C and Waters DD: Platelet hyperaggregability in patients with unstable angina. *Clin. Invest, Med II*: D77 1988.
16. Russell-Smith NC, Flower RF and Cardinal DC: Measuring platelet and leukocyte aggregation/adhesion responses in very small volumes of whole blood. *Jnl of Pharmac Meth*, 6, (315-333) 1981.
17. Selvaraj RJ, Sbarra AJ, Thomas GB, Cetrulo CL and Mitchell Jr GW: A microtechnique for studying chemiluminescence response of phagocytes using whole blood and its application to the evaluation of phagocytes in pregnancy. *RES: Jnl of the Reticuloendothelial Soc.* Vol. 31, (3-16) 1982
18. Ingerman-Wojenski et al
19. Greselle P, Zoja C, Deckmyn H, Arnout J, Vermynen J and Verstraete M: Dipyridamole Inhibits Platelet Aggregation in Whole Blood. *Thromb Haemostas* 50 (4) (852-856) 1983.
20. Riess et al
21. Russell-Smith et al
22. Challen A, Branch WJ and Cummings JH: Quantitation of platelet mass during aggregation in the electronic (Wellcome) whole blood aggregometer. *Jnl of Pharmac Meth*, 8(115-122) 1982.
23. Del Maschio A, Evangelista A, Rajtar G, Chen AM, Cerletti C and de Gaetano G: Platelet activation by polymorphonuclear leukocytes exposed to chemo-tactic agents. *American Physiological Society*, 0363-6135/90, H870-H879, 1990.
24. Ibid
25. Lacoste et al
26. Latour et al
27. Cardinal DC and Flower RJ: The electronic aggregometer: a novel device for assessing platelet behavior in blood. *Jnl. of Pharmac. Meth.* 3 (135-158) 1980.
28. Riess et al
29. Nichols WL, Kaese SE, Gastineau DA, Otteman LA and Bowie EJW: Bernard-Soulier Syndrome: whole blood diagnostic assays of platelets. *Mayo Clin Proc.* 64 (522-530) 1989.
30. Ray MJ, Hawson GAT, Just SJE, McLachlan G. and O'Brien M. Relationship of Platelet Aggregation to Bleeding After Cardiopulmonary Bypass. *Ann Thorac Surg*, 1994;57:981-6.
31. Ray MJ, Marsh NA, Just SJE, Perrin EJ, O'Brien MF and Hawson GAT. Preoperative Platelet Dysfunction Increases the Benefit of Aprotinin in Cardiopulmonary Bypass. *Ann Thorac Surg*, 1997;63:57-63
32. Mascelli MA, Worley S, Veriabo NJ, Lance ET, Mack S, Schaible X, Weisman HF and Jordan RE. Rapid Assessment of Platelet Function With a Modified Whole-Blood Aggregometer in Percutaneous Transluminal Coronary Angioplasty Patients Receiving Anti-GPIIb/IIIa Therapy. *Circulation* Vol 96, No. 11 December 2, 1997.
33. Stewart MW, Eches WX, Boshkov LK and Pordon PA. Heparin-Induced Thrombocytopenia: an Improved Method of Detection Based on Lumi-Aggregometry. *Brit Jnl of Haem*, 1995, 91, 173-177.
34. Bednar MM, Dooley RH, Tapanes R, Lublin JC and Gross CE: Ticlopidine Augments Lumino-Dependent Chemiluminescence in Human Neutrophils. *J Biol Chem* 1995.
35. Bednar MM, Dooley RH, Zamani M, Howard DB, Gross CE: Neutrophil and Platelet Activity and Quantification Following Delayed tPA Therapy in a Rabbit Model of Thromboembolic Stroke. *Journal of Thrombosis and Thrombolysis* 1995;1 (179-185)
36. Sathiropas P, Marbet GA, Sahaphong S, Duckert F: Detection of Small Inhibitory Effects of Acetylsalicylic Acid (ASA) by Platelet Impedance Aggregometry in Whole Blood. *Thrombosis Research* 51 (55-62) 1988.
37. Manoharan A, Gemmel R, Brighton T, Dunkley S, Lopez K and Kyle R: Thrombosis and bleeding in myeloproliferative disorders: identification of at-risk patients with whole blood platelet aggregation studies. *Brit Jnl of Haematology*, 1999, 105, 618-625.
38. Friend M, Vucenik I, Miller M: Reduced Platelet Responsiveness to Aspirin in Hyperlipidemic Subjects: Is 325 mg/d Sufficient? Poster Presented at American Heart Association 73rd Scientific Sessions [Abstract]
39. Friend M, Vucenik I, Miller M: Platelet responsiveness to aspirin in patients with hyperlipidaemia. *BMJ*, 2003;326: 83-83
40. Gu J, Prastein D, Pierson RN, White C, Griffith BP, Gurbel P, Manchio J, Poston R: Aspirin Resistance: An Underrecognized Risk Factor in Early Vein Graft Thrombosis After Off-Pump Coronary Artery Bypass (OPCAB). Presented at American Heart Association, Nov 7-20, 2004
41. Poston R, Gu J, Manchio J, Lee A, Brown J, Gammie J, White C, Griffith BP: Platelet function tests predict bleeding and thrombotic events after off-pump coronary bypass grafting. *European Journal of Cardio-thoracic Surgery* 27 (2005) 584-591
42. Dyszkiewicz-Korpanty AM, Frenkel EP, Sarode R: Approach to the Assessment of Platelet Function: Comparison between Optical-based Platelet-rich Plasma and Impedance-based Whole Blood Platelet Aggregation Methods. *Clin Appl Thrombosis/Hemostasis* 11(1):25-35, 2005
43. Ashraf T, Ahmed M, Talpur S, Kundi A, Farugui AMA, Jaffery AH, Fareed A: Competency Profile of Locally Manufactured Clopidogrel Lowplatelet and Foreign Manufactured Clopidogrel Plavix in Patients of Suspected Ischemic Heart Disease (CLAP-IHD). *Journal of Pakistan Med Assoc* Vol:55, No 10, Oct 2005
44. Miller, Jonathan L: Platelet-Type von Willebrand's Disease. *Thrombosis and Hemostasis*, Vol 7, No. 4, 1985

WHEECON INSTRUMENTS PRIVATE LIMITED
NO,7, SEVENTH STREET,
GOPALAPURAM, CHENNAI 600 086
PHONE :044 2493 9274 2811 2996
FAX:044 2811 9553
EMAIL:wheecon@vsnl.com Web:www.wheecon.com

CHRONO-LOG

C O R P O R A T I O N

2 WEST PARK ROAD, HAVERTOWN, PA 19083-4691, USA

1-800-247-6665 • In PA 610-853-1130 • Fax 610-853-3972

e-mail chronolog@chronolog.com • www.chronolog.com